PHACO-TIPS AND I/A-INSTRUMENTS (MULTIPLE USE)

ENGLISH

PHACO-TIPS AND I/A-INSTRUMENTS (MULTIPLE USE)

INSTRUCTIONS FOR USE

These instructions for use are to be read and used only in conjunction with the manufacturer's instruction manual. Failure to consult the manufacturer's manual could result in improper use of this product.

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Put pressure on the tip and screw the phaco tip into the handpiece up to the stop position. Attention: Overtighten the tip can damage the multiple use / single use wrench as well as it can cause residue from the wrench on the surface of the tip.

TREATMENT ACCORDING TO EN ISO 17664

General principles: Before every use, all instruments must be cleaned, disinfected and sterilized; this applies particularly to the first use after delivery since all instruments are delivered unsterilized (cleaning and disinfection after removal of protective shipping packaging; sterilization). An effective cleaning and disinfection is an absolute necessity for effective sterilization. As part of your responsibility for the sterility of instruments during use, please ensure that only appropriate device- and product-specific procedures for cleaning/disinfection are used, that the devices used (disinfectant, sterilizer) are regularly maintained and inspected and that the va-

lided parameters are adhered to with every cycle. Only qualified personnel should perform the treatment. In addition, please observe the applicable legal regulations in your country as well as the hygiene regulations of the doctor's practice or hospital. This is especially true for the different guidelines regarding an effective inactivation of prions. **CLEANING AND DISINFECTION** **Principles according to EN ISO 15883:** For the cleaning and disinfection of instruments, a mechanical process (disinfectant) should be used, if possible. Due to its clearly reduced effectiveness, a manual procedure – even when using an ultrasound bath – should only be used if a mechanical process is not available. The use of a manual cleaning and disinfection procedure must be guaranteed by an additional product- and procedure-specific validation and is the user's responsibility. In both cases, pretreatment must be performed.

PRETREATMENT


Immediately after use (within no more than 1 hour), rough impurities must be removed from the products. To do so, use running water or a disinfectant solution; the disinfectant should be free of aldehyde (otherwise, possible fixing of contaminants), have a tested effectiveness (e.g., German Hygiene and Microbiology Society [DGHM], or FDA approval or CE marking), • a tested program for thermal disinfection (at least 10 minutes at 93 °C) is used (with chemical disinfection there is the danger of disinfectant residue on the instruments), if possible, • the program used is suitable for the instruments and has sufficient rinsing cycles, • for washing down, only sterile or low-germ (max. 10 germs/ml) and low-endotoxin (max. 0.25 endotoxin units /ml) water (e.g. Aqua purificata) is used (EN 285, Appendix B), • the air used for drying is filtered and • the disinfectant is regularly maintained and inspected. **In selecting the cleaning agent used, make sure that** • it is basically suitable for cleaning the instruments (e.g. contains no abrasive materials), • if no thermal disinfection is used, a suitable disinfectant with a tested effectiveness (e.g. DGHM, or FDA approval

MACHINE CLEANING/DISINFECTION

In selecting the disinfectant, make sure that • the disinfectant basically has a tested effectiveness (e.g. DGHM, or FDA approval or CE marking), • a tested program for thermal disinfection (at least 10 minutes at 93 °C) is used (with chemical disinfection there is the danger of disinfectant residue on the instruments), if possible, • the program used is suitable for the instruments and has sufficient rinsing cycles, • for washing down, only sterile or low-germ (max. 10 germs/ml) and low-endotoxin (max. 0.25 endotoxin units /ml) water (e.g. Aqua purificata) is used (EN 285, Appendix B), • the air used for drying is filtered and • the disinfectant is regularly maintained and inspected. **In selecting the cleaning agent used, make sure that** • it is basically suitable for cleaning the instruments (e.g. contains no abrasive materials), • if no thermal disinfection is used, a suitable disinfectant with a tested effectiveness (e.g. DGHM, or FDA approval

or CE marking) is also used and that it is compatible with the cleaning agent used and • the chemicals used are compatible with the instruments (see chapter "Material Durability"). The concentrations indicated by the manufacturer of the cleaning agent and, if necessary, of the disinfectant, must absolutely be complied with. **Sequence:** 1. Place the instruments in the disinfectant. While doing so, make sure that the instruments do not touch each other. 2. Start the program. 3. After the end of the program, remove the instruments from the disinfectant. 4. Inspect and pack the instruments as quickly as possible after removal (see "Inspection", "maintenance" and "packing" sections, if necessary after additional drying in a clean place).

Cleaning and disinfection device or other measures must ensure sufficient drying. **INSPECTION AND TESTING** Check all instruments after cleaning or cleaning/disinfection for corrosion, damaged surfaces, chipping and contamination and sift out damaged instruments. Instruments that are still dirty must be re-cleaned and disinfected. **MAINTENANCE** Instrument oil or other oil may not be used.

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MADE IN SWITZERLAND

Manufacturer:
MEDICEL AG
Dornierstrasse 11
CH-9423 Altenrhein
SWITZERLAND

Tel. +41 71 727 10 50
Fax +41 71 727 10 55
info@medicel.com
www.medicel.com

<p>be used; other sterilization procedures are not permitted.</p> <p>Steam sterilization: <ul style="list-style-type: none">fractioned vacuum process (with sufficient product drying). The use of the less effective gravitation process must be guaranteed by an additional product-, sterilizer- and process-specific validation, and is the responsibility of the user (if necessary, longer sterilization times are required). validated steam sterilizer pursuant to E DIN EN 13060 or DIN EN 285 <ul style="list-style-type: none">process validated according to DIN EN/ANSI AAMI ISO 17665 (valid commissioning and product-specific performance appraisal) Sterilization temperature 134 °C +4/-13 °C (274 F° +6/-24 F°); with saturated steam pressure of 3.41 bar. Sterilization time (exposure time at sterilization temperature) of at least 20 minutes at 121 °C (250 F°) or 5 minutes at 134 °C (274 F°). <p>The flash sterilization process is fundamentally not permitted.Also, do not use any hot air sterilization, radiation sterilization, formalde-</p></p>	<p>hyde or ethylene oxide sterilization, and also no plasma sterilization.</p> <p>STORAGE</p> <p>After sterilization, the instruments must be stored dry.</p> <p>MATERIAL DURABILITY</p> <p>When choosing the cleansing agents and disinfectants, make sure that they do not contain the following ingredients: <ul style="list-style-type: none">strongacids orlyes aluminum chloride <ul style="list-style-type: none">halogens halogenated hydrocarbons furfural <ul style="list-style-type: none">phenols methylene chloride nitrobenzene <p>Never clean any of the instruments, sterilization trays and sterilization containers with metal brushes or steel wool. The permitted cleansing agents may only be used individually. All instruments, sterilization trays and sterilization containers may not be exposed to temperatures higher than 138 °C (280 F°)!</p> <p>REUSABILITY</p> <p>The instruments can be reused as long as they are handled carefully and are undamaged and uncontaminated. The user takes</p></p>
---------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------	--------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------

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