ERGOJECT™

– top-loaded –

medicel

SWISS TECHNOLOGY FOR SURGERY

Instructions for use

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Mode d'emploi

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MADE IN SWITZERLAND

Manufacturer: **MEDICEL AG**
Dornierstrasse 11
CH-9423 Altenrhein
SWITZERLAND

Tel. +41 71 727 10 50
Fax +41 71 727 10 55
info@medicel.com
www.medicel.com

medicel

SWISS TECHNOLOGY FOR SURGERY

ENGLISH

Medicel ERGOJECT™ Single-use Lens Injection System for injection of ONE-PIECE FOLDABLE LENSES

DESCRIPTION
With the aid of the ERGOJECT™ (Fig. 1) injection systems, the majority of all foldable, 1-piece acrylic lenses with an optical diameter up to 6.5 millimetres can be injected, provided that the safe injection of the lens using the ERGOJECT™ injection system was validated in advance by the lens manufacturer as per ISO 11979. In contrast to tradi-

tionally, the cartridge of the ERGOJECT™ system is integrated into the injector. The cartridges of the ERGOJECT™ models are designed with a large loading chamber. This enables the lenses to be loaded quickly, precisely and easily from above.
APPLICATION
Folding and injection of 1-piece folding lenses in the capsular bag or the sulcus after extracapsular cataract extraction.
INSTRUCTIONS FOR USE
1. Open the blister in a sterile environment and remove the sterile ERGOJECT™ injector set.
2. Fill the cartridge nozzle with viscoelastic solution of low to moderate viscosity directly from the cartridge tip (Fig. 2a).
3. Now wet the bottom of the loading chamber with viscoelastic solution from the end of the loading chamber cover (Fig. 2b). Do not completely fill the chamber with viscoelastic solution, as this can move the lens during insertion. Also put one drop of viscoelastic solution on the silicone plunger tip.

4. When using certain hydrophobic lens materials, the cartridge tip and loading chamber should be completely rinsed out with a BSS afterwards in order to achieve maximum hydration. This step is generally not necessary when using hydrophilic lens materials. The exclusive use of a BSS solution or the use of viscoelastic solutions with high viscosity is explicitly not recommended.
5. The cartridges are designed with a large loading chamber. This enables the lenses to be loaded quickly and precisely. Hold the ERGOJECT™ system so that the flexible cartridge wing can be guided with the index finger of your left hand. Open the loading chamber sufficiently wide to be able to easily insert the lens (Fig. 3).
6. Position the lens in the middle of the loading chamber bet-

ween the silicone plunger tip and the cartridge wall. For lenses with c-loop haptics, both haptics should be placed against the optical system using the forceps. At the same time, press the wings of the loading chamber together slightly to keep the haptics in the position in which they were placed (Fig. 4).
7. Close the cartridge, at the same time exerting light pressure on the lens optical system with the forceps. This ensures that the lens folds in the right direction (Fig. 5). As soon as the "click-lock" mechanism clicks in place, the lens is safely loaded and ready for injection (Fig. 6).
8. Push the blue injector plunger forward until the rear push plate is flush against the injector housing (Fig. 7) or until the wheel of the ERGOJECT™ injector moves. If the lens has been loaded by the surgical staff, the injector can be passed to the surgeon in this position.
9. There is a small lever on the surface of the injector that allows you to choose whether to operate the injector in forward and backward mode (lever in middle position) or in forward only mode (lever turned to the left) (Fig. 8). In forward only mode, the backward movement of the plunger is disabled. This can be an advantage when the lens and silicone plunger tip generate a strong recoil that makes it more difficult to operate the impeller (especially with ERGOJECT™ models for small incisions). When forward only mode is activated, this is signalled by a faint click when the plunger is pushed forward. You can switch between the forward and backward mode and the forward only mode at any time during the injection process by simply adjusting the lever.
10. Insert the cartridge tip through the incision and push over the iris to the proximal edge of the pupil.
11. Using your index finger, pull the wheel of the ERGOJECT™ injector slowly in order to push the lens forward (Fig. 9). When forward only mode is activated, you can hear a faint click when the plunger is pushed forward.
12. Depending on the lens model, turning the injector slightly backward mode (lever in middle position) or in forward only

mode (lever turned to the left) (Fig. 8). In forward only mode, the backward movement of the plunger is disabled. This can be an advantage when the lens and silicone plunger tip generate a strong recoil that makes it more difficult to operate the impeller (especially with ERGOJECT™ models for small incisions). When forward only mode is activated, this is signalled by a faint click when the plunger is pushed forward. You can switch between the forward and backward mode and the forward only mode at any time during the injection process by simply adjusting the lever.
10. Insert the cartridge tip through the incision and push over the iris to the proximal edge of the pupil.
11. Using your index finger, pull the wheel of the ERGOJECT™ injector slowly in order to push the lens forward (Fig. 9). When forward only mode is activated, you can hear a faint click when the plunger is pushed forward.
12. Depending on the lens model, turning the injector slightly backward mode (lever in middle position) or in forward only

into the capsular bag and then withdraw the injector tip from the backward movement of the plunger is disabled. This can be an advantage when the lens and silicone plunger tip generate a strong recoil that makes it more difficult to operate the impeller (especially with ERGOJECT™ models for small incisions). When forward only mode is activated, this is signalled by a faint click when the plunger is pushed forward. You can switch between the forward and backward mode and the forward only mode at any time during the injection process by simply adjusting the lever.
10. Insert the cartridge tip through the incision and push over the iris to the proximal edge of the pupil.
11. Using your index finger, pull the wheel of the ERGOJECT™ injector slowly in order to push the lens forward (Fig. 9). When forward only mode is activated, you can hear a faint click when the plunger is pushed forward.
12. Depending on the lens model, turning the injector slightly backward mode (lever in middle position) or in forward only

into the capsular bag and then withdraw the injector tip from the backward movement of the plunger is disabled. This can be an advantage when the lens and silicone plunger tip generate a strong recoil that makes it more difficult to operate the impeller (especially with ERGOJECT™ models for small incisions). When forward only mode is activated, this is signalled by a faint click when the plunger is pushed forward. You can switch between the forward and backward mode and the forward only mode at any time during the injection process by simply adjusting the lever.
10. Insert the cartridge tip through the incision and push over the iris to the proximal edge of the pupil.
11. Using your index finger, pull the wheel of the ERGOJECT™ injector slowly in order to push the lens forward (Fig. 9). When forward only mode is activated, you can hear a faint click when the plunger is pushed forward.
12. Depending on the lens model, turning the injector slightly backward mode (lever in middle position) or in forward only

product defects, which are clearly not the result of incorrect handling or the use of lenses not validated with this injector model.
ATTENTION: US federal regulations restrict the sale of this product to medical practitioners and those acting on their behalf.*

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








ENGLISH

Medicel ERGOJECT™

Single-use Lens Injection System for injection of ONE-PIECE FOLDABLE LENSES

DESCRIPTION

With the aid of the ERGOJECT™ (Fig. 1) injection systems, the majority of all foldable, 1-piece hydrophobic acrylic lenses with an optical diameter up to 6.5 millimetres can be injected, provided that the safe injection of the lens using the respective ERGOJECT™ injection system was validated in advance by the lens manufacturer as per ISO 11979. This system is not

suitable for injecting hydrophilic lenses. In contrast to traditional injection systems, the cartridge of the ERGOJECT™ system is integrated into the injector. The cartridges of the ERGOJECT™ models are designed with a large covered loading chamber. This enables the lenses to be loaded quickly, precisely and easily. APPLICATION Folding and injection of 1-piece hydrophobic folding lenses in the capsular bag or the sulcus after extracapsular cataract extraction. INSTRUCTIONS FOR USE 1. Open the blister in a sterile environment and remove the sterile ERGOJECT™ injector set. 2. Fill the cartridge nozzle with viscoelastic solution of low to moderate viscosity directly from the cartridge tip (Fig. 2a). 3. Now wet the bottom of the loading chamber with viscoelastic solution from the end of the loading chamber cover (Fig. 2b). Do not completely fill the chamber with viscoelastic solution, as this can lead to the lens during insertion. Also put one drop of viscoelastic solution on the silicone plunger tip. 4. When using certain hydrophobic lens materials, the cartridge tip and loading chamber should be completely rinsed out with a BSS afterwards in order to achieve maximum hydration. This step is generally not necessary when using hydrophilic lens materials. The exclusive use of a BSS solution or the use of viscoelastic solutions with high viscosity is explicitly not recommended. The lens should be injected 3 minutes at the latest after filling with the viscoelastic solution. Viscoelastic substances can lose their lubrication properties if they are in contact with air for too long. 5. The cartridges are designed with a large loading chamber. This enables the lenses to be loaded quickly and precisely. Hold the ERGOJECT™ system so that the flexible cartridge wing can be guided with the index finger of your left hand (Fig. 3).

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