GÜELL DMEK

Instructions for use

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Mode d'emploi

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MADE IN SWITZERLAND

Manufacturer: **MEDICEL AG**
Dornierstrasse 11
9423 Altenrhein
SWITZERLAND

Tel. +41 71 727 10 50
Fax +41 71 727 10 55
info@medicel.com
www.medicel.com

ENGLISH

Medicel GÜELL DMEK
Single-use corneal endothelium injection system

for the injection of
ENDOTHELIAL IMPLANTS IN
DESCEMET MEMBRANE ENDOTHELIAL
KERATOPLASTY (DMEK)

DESCRIPTION
Thanks to the GÜELL DMEK hydraulic system, the GÜELL DMEK injection system allows implantation of corneal transplants using Descemet Membrane Endothelial Keratoplasty (DMEK).

The GÜELL DMEK injection system consists of the sterile GÜELL DMEK single-use injector with silicone tip.

APPLICATION
Injection of self-rolling corneal endothelium transplants in Descemet Membrane Endothelial Keratoplasty (DMEK). Non-self-rolling corneal endothelium implants with a high stroma proportion, such as those used in Descemet's Stripping Automated Endothelial Keratoplasty (DSAEK), are not suitable for implantation with the GÜELL DMEK injection system. The system is also unsuitable for the injection of intraocular lenses in cataract operations.

INSTRUCTIONS FOR USE

1. Open the blister in a sterile environment and remove the injector (Fig. 1).
2. Fill the cartridge tip and the loading chamber of the GÜELL DMEK injector with BSS or organic culture medium until the cartridge tip and loading chamber are completely filled and covered with BSS respectively organic culture medium (Fig. 2).
3. Gently lift the rolled implant from the tray, using sterile flat forceps and place in the loading chamber (Fig. 3). Minimize contact between the endothelium and the outer walls of the loading chamber. Next, flush the rolled implant along of the loading chamber into the cartridge tip using BSS, avoiding direct contact with the implant (Fig. 4), until the implant lies entirely within the cartridge tip. Depending on the surgical technique, an air bubble can be placed behind the implant in the cartridge tip (Fig. 5).
4. Then squeeze the loading chamber clicks together until the "click-lock" mechanism clicks into place.
5. Push the blue injector plunger forward until the rear push

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