Instructions for use Accuject 3P

| | | | |
|----|---------------------------|---------------------|----|
| EN | Instructions for use | Page | 2 |
| DE | Gebrauchsanweisung | Seite | 5 |
| AR | تعليمات الاستخدام | صفحة | 8 |
| BG | Инструкции за употреба | страница на книгата | 11 |
| CS | Pokyny k použití | S stránka knihy | 14 |
| DA | Brugsanvisning | Bog side | 17 |
| EL | Οδηγίες χρήσης | σελίδα βιβλίου | 20 |
| ES | Instrucciones de uso | Página | 23 |
| ET | Kasutusjuhend | Raamatu leht | 26 |
| FI | Käyttöohjeet | Kirjan sivu | 29 |
| FR | Mode d'emploi | Page | 32 |
| HR | Upute za upotrebu | Stranica knjige | 35 |
| HU | Használati utasítás | Könyv oldala | 38 |
| IT | Istruzioni per l'uso | Pagina | 41 |
| KO | 사용 지침 | 예약 페이지 | 44 |
| LT | Naudojimo instrukcija | Knygos puslapis | 47 |
| LV | Lietošanas instrukcija | Grāmatas lapa | 50 |
| NL | Gebruiksaanwijzing | Pagina | 53 |
| NO | Bruksanvisning | Bokside | 56 |
| PL | Instrukcja użytkowania | Strona książki | 59 |
| PT | Instruções de utilização | Página | 62 |
| RO | Instrucțiuni de utilizare | Pagina | 65 |
| RU | Инструкция по применению | Страница книги | 68 |
| SK | Návod na použitie | Stránka knihy | 71 |
| SL | Navodila za uporabo | Stran knjige | 74 |
| SR | Uputstvo za upotrebu | Stranica knjige | 77 |
| SV | Bruksanvisning | Sida | 80 |
| TR | Kullanma talimatları | Kitap sayfası | 83 |
| UK | Інструкція з експлуатації | Сторінка книги | 86 |

1. Scope of application

These instructions for use apply to the following medical devices:

| Reference number | Product name |
|------------------|------------------------------|
| LP604520 | Accuject 2.5-3P Injector Set |
| LP604580 | Accuject 3.0-3P Injector Set |

2. Scope

Intended use: Folding and injecting of an intraocular lens (IOL) during eye surgery.

Indication: Visual problems that require implantation of IOL.

Contraindications: There are no specific contraindications for the injectors, apart from general contraindications related to ophthalmic surgery.

  User Group: The medical device must only be handled by healthcare professionals and applied to the patient by the surgeon. Use by unqualified personnel could result in patient injury.

 Patient Group: A patient-related, careful preoperative assessment, with well-founded clinical judgment and risk-benefit assessment, is the responsibility of the surgeon.

Clinical benefit: IOL injector allows IOL implantation through small incisions.

3. Notes on safety

The suitability of the intraocular lens (IOL) in combination with the medical device (IOL injector) must be tested and approved in advance by the intraocular lens manufacturer.

  **Warning**
Using a damaged or non-sterile medical device can lead to infection / endophthalmitis in the patient.
Do not use the medical device if there are signs of damage to the medical device itself or to the sterile packaging.

 **Warning**
The medical device must not be used after the expiry date. Use after the expiration date may result in infection / endophthalmitis of the patient.

 **Warning**
Touching the endothelium with the silicone cushion or IOL can damage the corneal endothelium.

 **Warning**
Appropriate surgical technique is the responsibility of the individual surgeon. The surgeon has to assess the suitability of the respective procedure on the basis of his / her training and experience.

 **Warning**
The medical device is intended for single use only. Reuse or inappropriate reprocessing may lead to serious adverse effects on the health and safety of the patient.

  - Reuse may result in infection / endophthalmitis, TASS or allergic reaction due to cleaning residues.
- The performance of the medical device may be compromised. Heat or the cleaning process will destroy performance properties (e.g. mechanical properties and gliding properties) of the medical device.

 **Caution**
In order to avoid capsule rupture, the IOL must be injected slowly and carefully into the eye.

4. Reporting to manufacturer and authorities

Serious incidents occurring in connection with the device must be reported to the manufacturer and the competent authority in which the user and/or patient is established.

5. Application instructions for the single-use IOL injector


1. Open the blister in a sterile environment and remove the sterile injector (1).
Make sure that the injector plunger is not advanced prematurely.
2. Hold the injector so that the rear wing of the loading chamber can be guided with the index finger of your left hand (2).
3. Fill the cartridge tip and the loading chamber sufficiently with ophthalmic viscosurgical devices (OVD) (3).
It is recommended to flush the cartridge tip and the loading chamber with Balanced Salt Solution (BSS) before the application of OVD.

Allow the OVD to react (30s). OVD can lose their lubricating properties if they are in contact with air for too long. Therefore, the IOL should be injected without delay after loading.

Hydrophobic IOL types: The application of OVD is required. The exclusive use of BSS is not permitted.
4. Open the loading chamber so that the IOL can be easily loaded.

Position the lens in the middle of the loading chamber with a concave orientation and make sure that the rear haptic comes to rest over the green injector plunger head (4). Make sure that the front haptics are positioned under the positioning hook.
5. Use the sterile forceps to apply slight pressure to the IOL optics (5a), this ensures that the IOL folds in a concave direction.
6. At the same time, close the wings of the loading chamber until the „Click-Lock“ mechanism engages (5b).
7. Loosen the attachment on the front part of the injector and slowly pull the purple positioning hook out of the cartridge tip. This will automatically stretch the front haptics of the lens.
8. Advance the IOL into the conical tip of the cartridge (6).

The IOL is now loaded and ready for injection.
9. Push the injector plunger forward. As the distal haptic emerges, turn the injector slightly clockwise to allow the haptic to fold out into a flat profile.
10. With the distal haptic in the correct position, slowly inject the IOL and when the optic emerges from the cartridge tip, turn it counterclockwise back to its original position. Relieve the injector plunger until the proximal haptic is released from the plunger. In a second movement, implant the proximal haptic directly into the capsular bag.
11. If necessary, assist the IOL with the help of a suitable positioning hook during the exit and bring the IOL into its final position.
12. Thoroughly remove the viscoelastic material from the eye and IOL using standard irrigation and aspiration techniques.

6. Disposal

Warning

The medical device may be contaminated with potentially infectious substances of human origin after use.

Dispose of the medical device and packaging after use according to the applicable guidelines of biological hazardous waste.

7. Symbols

| | | | |
|--|---|--|---|
| | Sterilized using ethylene oxide | | Keep dry |
| | Number of medical device in packaging unit | | Do not use if package is damaged and consult instructions for use |
| | Single sterile barrier system | | Keep away from sunlight |
| | Single sterile barrier system with protecting packaging outside | | Do not re-sterilize |
| | Date of manufacture YYYY-MM-DD and country of manufacture | | CE symbol with number of the notified body |
| | Do not re-use | | Authorized representative in the European Community / European Union |
| | Reference number | | Open here |
| | Lot number | | Unique device identifier |
| | Expiry date | | Medical device |
| | Manufacturer | | Consult electronic instructions for use www.medicel.com/ifu |
| | To be used by qualified personnel only | | Caution |
| | Recycling | | |

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