## Instructions for use Accuject Dual Back-Loaded

EN	Instructions for use	Page	2
DE	Gebrauchsanweisung	Seite	5
AR	إدخال تعليمات	باتة ءءءفء	8
BG	Инструкции за употреба	страница на книгата	11
CS	Pokyny k použití	Sstránka knihy	14
DA	Brugsanvisning	Bog side	17
EL	Οδηγίες χρήσης	σελίδα βιβλίου	20
ES	Instrucciones de uso	Página	23
ET	Kasutusjuhend	Raamatu leht	26
FI	Käyttöohjeet	Kirjan sivu	29
FR	Mode d'emploi	Page	32
HR	Upute za upotrebu	Stranica knjige	35
HU	Használati utasítás	Könyv oldala	38
IT	Istruzioni per l'uso	Pagina	41
KO	사용 지침	예약 페이지	44
LT	Naudojimo instrukcija	Knygos puslapis	47
LV	Lietošanas instrukcija	Grāmatas lapa	50
NL	Gebruiksaanwijzing	Pagina	53
NO	Bruksanvisning	Bokside	56
PL	Instrukcja użytkowania	Strona książki	59
PT	Instruções de utilização	Página	62
RO	Instrucţiuni de utilizare	Pagina	65
RU	Инструкция по применению	Страница книги	68
SK	Návod na použitie	Stránka knihy	71
SL	Navodila za uporabo	Stran knjige	74
SR	Uputstvo za upotrebu	Stranica knjige	77
SV	Bruksanvisning	Sida	80
TR	Kullanma talimatları	Kitap sayfası	83
UK	Інструкція з експлуатації	Сторінка книги	86

## 1. Scope of application

These instructions for use apply to the following medical devices:

Reference number	Product name
ADB1600	Accuject Dual 1.6-BL Injector
ADB1800	Accuject Dual 1.8-BL Injector
ADB2000	Accuject Dual 2.0-BL Injector
ADB2100	Accuject Dual 2.1-BL Injector
ADB2200	Accuject Dual 2.2-BL Injector
ADB2400	Accuject Dual 2.4-BL Injector
ADB2600	Accuject Dual 2.6-BL Injector
ADB3000	Accuject Dual 3.0-BL Injector

## 2. Intended field of application

The suitability of the intraocular lens (IOL) in combination with the medical device (IOL injector) must be tested and approved in advance by the intraocular lens manufacturer.

### 2.1 Intended use

Folding and injecting of an intraocular lens (IOL) during eye surgery.

### 2.2 Clinical benefit

IOL injector allows IOL implantation through small incisions.

## 3. Notes on safety

### Warning

Using a damaged or non-sterile medical device can lead to infection / endophthalmitis in the patient. Do not use the medical device if there are signs of damage to the medical device itself or to the sterile packaging.

### Warning

The medical device must not be used after the expiry date. Use after the expiration date may result in infection / endophthalmitis of the patient.

### Warning

Touching the endothelium with the silicone cushion or IOL can damage the corneal endothelium.

### Warning

Appropriate surgical technique is the responsibility of the individual surgeon. The surgeon has to assess the suitability of the respective procedure on the basis of his / her training and experience.

### Warning

The medical device is intended for single use only. Reuse or inappropriate reprocessing may lead to serious adverse effects on the health and safety of the patient.

- Reuse may result in infection / endophthalmitis, TASS or allergic reaction due to cleaning residues.

- The performance of the medical device may be compromised. Heat or the cleaning process will destroy performance properties (e.g. mechanical properties and gliding properties) of the medical device.

### Warning

The medical device may be handled by healthcare professionals and applied to the patient by the surgeon. Use by unqualified personnel could result in patient injury.

A patient-related, careful preoperative assessment, with well-founded clinical judgment and risk-benefit assessment, is the responsibility of the surgeon. U.S. federal law restricts the sale of this medical device to physicians or on the order of physicians.

### Caution

In order to avoid capsule rupture, the IOL must be injected slowly and carefully into the eye.

#### 4. Reporting to manufacturer and authorities

Serious incidents occurring in connection with the device must be reported to the manufacturer and the competent authority in which the user and/or patient is established.

#### 5. Application instructions for the single-use IOL injector

1. Open the blister in a sterile environment and remove the sterile injector (1).
2. Hold the injector so that the rear wing of the loading chamber can be guided with the index finger of your left hand (2).
3. Fill the cartridge tip and, from the end of the loading chamber cover, the bottom of the loading chamber sufficiently with ophthalmic viscosurgical devices (OVD) (3).

Do not completely fill the loading chamber with OVD, as this can shift the IOL when it is inserted.

It is recommended to flush the cartridge tip and the loading chamber with Balanced Salt Solution (BSS) before the application of OVD.

Allow the OVD to react (30s). OVD can lose their lubricating properties if they are in contact with air for too long. Therefore, the IOL should be injected without delay after loading.

Hydrophobic IOL types: The application of OVD is required. The exclusive use of BSS is not permitted.

In addition, apply a drop of OVD to the silicone cushion (3).

The premature (even partial) compression of the loading chamber wings before positioning and loading the IOL can damage the injector and impair the functionality of the injector.

4. Position the IOL on the rear, uncovered part of the loading chamber (4). Press against the rear end of the optic with the sterile forceps. For IOLs with c-loop haptics, only press against the optics and not against the rear haptics.

Push the IOL into the loading chamber until the forceps touch the end of the guide slot in the cover (5).

5. Hold the forceps in this position and close the wings of the loading chamber with your second hand (6).

When using IOLs with c-loop haptics, do not pull out the forceps until the wings of the loading chamber are already slightly compressed and pressure is already being exerted on the IOL. This prevents the IOL from slipping back independently from its advanced position and the already pre-folded front haptic from stretching back to its original shape.

6. Close the wings of the loading chamber until the „Click-Lock“ mechanism engages (7).

7. Push the injector plunger forward with the finger flanges folded in (8).

If the silicone cushion is inclined and cannot be inserted into the closed loading chamber, the silicone cushion can be aligned with the forceps.

8. Advance the IOL into the conical tip of the cartridge as far as it will go or until half of the silicone cushion is visible in the viewing window (9).

Pull back the injector plunger a few millimetres and then push it forward again. This step ensures that the IOL is always gripped correctly.

The IOL is now loaded and ready for injection.


9. Syringe injector: unfold the finger flanges.  
Screw injector: leave the finger flange folded in.  
The finger flanges must be fully retracted for use in screw mode.

10. Slowly inject the IOL into the eye by either applying even pressure to the injector plunger (syringe injector) or slowly turning the injector plunger knob clockwise (screw injector) (10).

11. The correct exit of the IOL can be supported by turning the injector slightly. Only push the injector plunger forward until the IOL has emerged completely.

12. If necessary, assist the IOL with the help of a suitable positioning hook during the exit and bring the IOL into its final position.

13. Thoroughly remove the viscoelastic material from the eye and IOL using standard irrigation and aspiration techniques.

#### 6. Disposal

##### Warning

The medical device may be contaminated with potentially infectious substances of human origin after use. Dispose of the medical device and packaging after use according to the applicable guidelines of biological hazardous waste.

#### 7. Symbols

	Sterilized using ethylene oxide		Keep dry
	Number of medical device in packaging unit		Do not use if package is damaged and consult instructions for use
	Single sterile barrier system		Keep away from sunlight
	Single sterile barrier system with protecting packaging outside		Do not re-sterilize
	Date of manufacture YYYY-MM-DD and country of manufacture		CE symbol with number of the notified body
	Do not re-use		Authorized representative in the European Community / European Union
	Reference number		Open here
	Lot number		Unique device identifier
	Expiry date		Medical device
	Manufacturer		Consult electronic instructions for us <a href="http://www.medicel.com/ifu">www.medicel.com/ifu</a>
	To be used by qualified personnel only		Caution
	Recycling		

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Referenznummer	Name des Produkts
ADB1600	Accuject Dual 1.6-BL Injector
ADB1800	Accuject Dual 1.8-BL Injector
ADB2000	Accuject Dual 2.0-BL Injector
ADB2100	Accuject Dual 2.1-BL Injector
ADB2200	Accuject Dual 2.2-BL Injector
ADB2400	Accuject Dual 2.4-BL Injector
ADB2600	Accuject Dual 2.6-BL Injector
ADB3000	Accuject Dual 3.0-BL Injector

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Referentni broj	Ime proizvoda
ADB1600	Accuject Dual 1.6-BL Injector
ADB1800	Accuject Dual 1.8-BL Injector
ADB2000	Accuject Dual 2.0-BL Injector
ADB2100	Accuject Dual 2.1-BL Injector
ADB2200	Accuject Dual 2.2-BL Injector
ADB2400	Accuject Dual 2.4-BL Injector
ADB2600	Accuject Dual 2.6-BL Injector
ADB3000	Accuject Dual 3.0-BL Injector

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
Referensnummer	Produktnamn
ADB1600	Accuject Dual 1.6-BL Injector
ADB1800	Accuject Dual 1.8-BL Injector
ADB2000	Accuject Dual 2.0-BL Injector
ADB2100	Accuject Dual 2.1-BL Injector
ADB2200	Accuject Dual 2.2-BL Injector
ADB2400	Accuject Dual 2.4-BL Injector
ADB2600	Accuject Dual 2.6-BL Injector
ADB3000	Accuject Dual 3.0-BL Injector

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