Instructions for use Ergoject PRO

| | | | |
|----|---------------------------|---------------------|----|
| EN | Instructions for use | Page | 2 |
| DE | Gebrauchsanweisung | Seite | 5 |
| AR | تعليمات الاستخدام | صفحة | 8 |
| BG | Инструкции за употреба | страница на книгата | 11 |
| CS | Pokyny k použití | S stránka knihy | 14 |
| DA | Brugsanvisning | Bog side | 17 |
| EL | Οδηγίες χρήσης | σελίδα βιβλίου | 20 |
| ES | Instrucciones de uso | Página | 23 |
| ET | Kasutusjuhend | Raamatu leht | 26 |
| FI | Käyttöohjeet | Kirjan sivu | 29 |
| FR | Mode d'emploi | Page | 32 |
| HR | Upute za upotrebu | Stranica knjige | 35 |
| HU | Használati utasítás | Könyv oldala | 38 |
| IT | Istruzioni per l'uso | Pagina | 41 |
| KO | 사용 지침 | 예약 페이지 | 44 |
| LT | Naudojimo instrukcija | Knygos puslapis | 47 |
| LV | Lietošanas instrukcija | Grāmatas lapa | 50 |
| NL | Gebruiksaanwijzing | Pagina | 53 |
| NO | Bruksanvisning | Bokside | 56 |
| PL | Instrukcja użytkowania | Strona książki | 59 |
| PT | Instruções de utilização | Página | 62 |
| RO | Instrucţiuni de utilizare | Pagina | 65 |
| RU | Инструкция по применению | Страница книги | 68 |
| SK | Návod na použitie | Stránka knihy | 71 |
| SL | Navodila za uporabo | Stran knjige | 74 |
| SR | Uputstvo za upotrebu | Stranica knjige | 77 |
| SV | Bruksanvisning | Sida | 80 |
| TR | Kullanma talimatları | Kitap sayfası | 83 |
| UK | Інструкція з експлуатації | Сторінка книги | 86 |

1. Scope of application

These instructions for use apply to the following medical devices:

| Reference number | Product name |
|------------------|------------------------------|
| ERP1600 | Ergoject PRO 1.6-TL Injector |
| ERP1800 | Ergoject PRO 1.8-TL Injector |
| ERP2000 | Ergoject PRO 2.0-TL Injector |
| ERP2100 | Ergoject PRO 2.1-TL Injector |
| ERP2200 | Ergoject PRO 2.2-TL Injector |
| ERP2400 | Ergoject PRO 2.4-TL Injector |
| ERP2600 | Ergoject PRO 2.6-TL Injector |
| ERP3000 | Ergoject PRO 3.0-TL Injector |

2. Intended field of application

The suitability of the intraocular lens (IOL) in combination with the medical device (IOL injector) must be tested and approved in advance by the intraocular lens manufacturer.

2.1 Intended use

Folding and injecting of an intraocular lens (IOL) during eye surgery.

2.2 Clinical benefit

IOL injector allows IOL implantation through small incisions.

3. Notes on safety

Warning

Using a damaged or non-sterile medical device can lead to infection / endophthalmitis in the patient. Do not use the medical device if there are signs of damage to the medical device itself or to the sterile packaging.

Warning

The medical device must not be used after the expiry date. Use after the expiration date may result in infection / endophthalmitis of the patient.

Warning

Touching the endothelium with the silicone cushion or IOL can damage the corneal endothelium.

Warning

Appropriate surgical technique is the responsibility of the individual surgeon. The surgeon has to assess the suitability of the respective procedure on the basis of his / her training and experience.

Warning

The medical device is intended for single use only. Reuse or inappropriate reprocessing may lead to serious adverse effects on the health and safety of the patient.

- Reuse may result in infection / endophthalmitis, TASS or allergic reaction due to cleaning residues.

- The performance of the medical device may be compromised. Heat or the cleaning process will destroy performance properties (e.g. mechanical properties and gliding properties) of the medical device.

Warning

The medical device may be handled by healthcare professionals and applied to the patient by the surgeon. Use by unqualified personnel could result in patient injury.

A patient-related, careful preoperative assessment, with well-founded clinical judgment and risk-benefit assessment, is the responsibility of the surgeon. U.S. federal law restricts the sale of this medical device to physicians or on the order of physicians.

Caution

In order to avoid capsule rupture, the IOL must be injected slowly and carefully into the eye.

4. Reporting to manufacturer and authorities

Serious incidents occurring in connection with the device must be reported to the manufacturer and the competent authority in which the user and/or patient is established.

5. Application instructions for the single-use IOL injector

1. Open the blister in a sterile environment and remove the sterile injector (1).
2. Open the aluminium cap of the lens case in a sterile environment.
3. Remove the loading chamber with the lens from the lens case by pulling up on the upper loop of the lens holder (2).
4. Insert the loading chamber with the lens into the opening in the injector body and press it in until it audibly clicks into place (3).
5. Fill the cartridge tip and the bottom of the loading chamber sufficiently with ophthalmic viscosurgical devices (OVD) (4).

Do not completely fill the loading chamber with OVD, as this can shift the IOL when it is inserted.

It is recommended to flush the cartridge tip and the loading chamber with Balanced Salt Solution (BSS) before the application of OVD.

Allow the OVD to react (30s). OVD can lose their lubricating properties if they are in contact with air for too long. Therefore, the IOL should be injected without delay after loading.

The application of a drop of OVD onto the silicone cushion is required.

6. Hold the loop of the lens holder with your index and middle finger and use your thumb to gently flip up the lower hinge of the lens holder (5). This unlatches the hinge and allows the lens holder to be pulled up and removed without any force.

Do not pull the lens holder away to the side, but only remove it vertically upwards, since pulling the lens holder sideways can shift the lens in the loading chamber. Hold the injector horizontally during and after removal of the lens holder to ensure the lens does not slip within the loading chamber.

7. close the wings of the loading chamber until the „Click-Lock“ mechanism engages (6).
8. Push the injector plunger forward.

If the silicone cushion is inclined and cannot be inserted into the closed loading chamber, the silicone cushion can be aligned with the forceps.

Pull back the injector plunger a few millimetres and then push it forward again. This step ensures that the IOL is always gripped correctly.


9. Push the injector plunger forward until the rear pusher plate is flush with the injector housing or until the wheel moves (7).
The IOL is now loaded and ready for injection.

10. Use the small lever on the injector housing to choose whether the injector should be operated in forward and reverse mode (lever in middle position) or in forward mode only (lever to the left) (8). In the forward mode only, the return movement of the injector plunger is blocked. This is advantageous if the IOL and silicone cushion generate a strong recoil, which makes it difficult to operate the injector wheel (especially with models for small incisions). The activation of the exclusive advance mode is signalled by a slight click when the injector plunger is pushed forward. At any time during the injection process, it is possible to switch back and forth between forward and reverse mode and exclusive forward mode by simply adjusting the lever.

11. With your index finger, slowly pull the wheel back to advance the IOL (9).

If the rear haptic is wedged between the injector plunger and the cartridge tip after the plunger has been fully advanced, set the injector to forward and reverse mode (lever in the middle position) and pull the injector plunger back a few millimetres until the haptic is exposed. Activate exclusive forward mode again if necessary and inject the rear haptics in a second push.

12. Slowly inject the IOL into the eye.
13. The correct exit of the IOL can be supported by turning the injector slightly. Only push the injector plunger forward until the IOL has emerged completely.
14. If necessary, assist the IOL with the help of a suitable positioning hook during the exit and bring the IOL into its final position.
15. Thoroughly remove the viscoelastic material from the eye and IOL using standard irrigation and aspiration techniques.

6. Disposal

Warning

The medical device may be contaminated with potentially infectious substances of human origin after use.
Dispose of the medical device and packaging after use according to the applicable guidelines of biological hazardous waste.

7. Symbols

| | | | |
|--|---|--|--|
| | Sterilized using ethylene oxide | | Keep dry |
| | Number of medical device in packaging unit | | Do not use if package is damaged and consult instructions for use |
| | Single sterile barrier system | | Keep away from sunlight |
| | Single sterile barrier system with protecting packaging outside | | Do not re-sterilize |
| | Date of manufacture YYYY-MM-DD and country of manufacture | | CE symbol with number of the notified body |
| | Do not re-use | | Authorized representative in the European Community / European Union |
| | Reference number | | Open here |
| | Lot number | | Unique device identifier |
| | Expiry date | | Medical device |
| | Manufacturer | | Consult electronic instructions for us www.medicel.com/ifu |
| | To be used by qualified personnel only | | Caution |
| | Recycling | | |

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