I/A Handpieces Coaxial and Bi-Manual

Instructions for use

 Medicel AG
Dornierstrasse 11

9423 Altenrhein
Switzerland

t +41 71 727 10 50
f +41 71 727 10 55

info@medicel.com
www.medicel.com

1

2

3

4

5

EN	Instructions for use	2
DE	Gebrauchsanweisung	4
AR	تعليمات الاستخدام	6
BG	Инструкции за употреба	8
CS	Pokyny k použití	10
DA	Brugsanvisning	12
EL	Οδηγίες χρήσης	14
ES	Instrucciones de uso	16
ET	Kasutusjuhend	18
FI	Käyttöohjeet	20

FR	Mode d'emploi	22
HR	Upute za upotrebu	24
HU	Használati utasítás	26
IT	Istruzioni per l'uso	28
KO	사용 지침	30
LT	Naudojimo instrukcija	32
LV	Lietošanas instrukcija	34
NL	Gebruiksaanwijzing	36
NO	Bruksanvisning	38
PL	Instrukcja użytkowania	40

PT	Instruções de utilização	42
RO	Instrucțiuni de utilizare	44
RU	Инструкция по применению	46
SK	Návod na použitie	48
SL	Navodila za uporabo	50
SR	Uputstvo za upotrebu	52
SV	Bruksanvisning	54
TR	Kullanma talimatları	56
UK	Інструкція з експлуатації	58

1. Scope of application

Instructions for use apply to: • BICEL terralution • Swiss Bi-Manual I/A System
• Bi-Manual I/A • Swiss Coaxial I/A System • Coaxial I/A
• Femto I/A Handpieces

2. Scope

Intended use: Removal of human lens during eye surgery.

Indication: Visual problems that require removal of lens.

Contraindications: There are no specific contraindications for the medical device, apart from general contraindications related to ophthalmic surgery.

User group: The medical device must only be handled by healthcare professionals and applied to the patient by the surgeon.

Use by unqualified personnel could result in patient injury.

Patient group: A patient-related, careful preoperative assessment, with well-founded clinical judgment and risk-benefit assessment, is the responsibility of the surgeon.

Clinical benefit: Removal of human lens during eye surgery through small incisions.

Rx ONLY

3. Notes on safety

Warning

Using a damaged or non-sterile medical device can lead to an infection or endophthalmitis of the patient. Do not use if package is damaged or sterile barrier is harmed.

Warning

Avoid contact with sharp and / or pointed device components in order to facilitate safe use and disposal of the device.

Warning

The medical device must not be used after its expiry date. Use after the expiration date may result in infection / endophthalmitis of the patient.

Warning

Never reuse the medical device. Reuse and / or reesterilization may compromise the device performance and safety, which could cause infections and serious harm to the patient's health and safety.

Never reprocess and re-use the device.

Warning

The medical device may be handled by healthcare professionals and applied to the patient by the surgeon. Use by unqualified personnel could result in patient injury. A patient-related, careful preoperative assessment, with well-founded clinical judgment and

Rx ONLY

risk-benefit assessment, is the responsibility of the surgeon.

Caution

Make sure to always combine Mediceal components of the correct size.

Caution

The modification on the instruments e.g. bending of the tip is not allowed as this could lead to instrument failure (break)

Caution

The Irrigation and Aspiration connectors need to be attached to the handpiece with a twist of the connectors.

4. Application instructions for the single-use I/A Handpieces

Caution

See illustrated instruction for use and attach the Irrigation and Aspiration Connectors with a twist of the connectors. Read and observe the respective surgical system's instructions for use.

5. Disposal

Warning

The medical device may be contaminated with potentially infectious substances of human origin after use. Dispose of the medical device and packaging after use according to the applicable guidelines of biological hazardous waste.

6. Notifications to manufacturer and authorities

Serious incidents occurring in connection with the device must be reported to the manufacturer and the competent authority in which the user and /or patient is established.

7. Symbols not defined by ISO 15223-1:

	Number of medical device in packaging unit		Open here
	To be used by qualified personnel only		Recycling

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Gebrauchsanweisung gilt für: • BICEL terralution • Swiss Bi-Manual I/A System
• Bi-Manual I/A • Swiss Coaxial I/A System • Coaxial I/A
• Femto I/A Handpieces

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