Phaco Accessories Sets single-use

Instructions for use


 Medicel AG
Dornierstrasse 11

9423 Altenrhein
Switzerland

t +41 71 727 10 50
f +41 71 727 10 55

info@medicel.com
www.medicel.com

1

2

3

4

5

6

7

8

EN	Instructions for use	2	FR	Mode d'emploi	22	PT	Instruções de utilização	42
DE	Gebrauchsanweisung	4	HR	Upute za upotrebu	24	RO	Instrucțiuni de utilizare	44
AR	تعليمات الاستخدام	6	HU	Használati utasítás	26	RU	Инструкция по применению	46
BG	Инструкции за употреба	8	IT	Istruzioni per l'uso	28	SK	Návod na použitie	48
CS	Pokyny k použití	10	KO	사용 지침	30	SL	Navodila za uporabo	50
DA	Brugsanvisning	12	LT	Naudojimo instrukcija	32	SR	Uputstvo za upotrebu	52
EL	Οδηγίες χρήσης	14	LV	Lietošanas instrukcija	34	SV	Bruksanvisning	54
ES	Instrucciones de uso	16	NL	Gebruiksaanwijzing	36	TR	Kullanma talimatları	56
ET	Kasutusjuhend	18	NO	Bruksanvisning	38	UK	Інструкція з експлуатації	58
FI	Käyttöohjeet	20	PL	Instrukcja użytkowania	40			

1. Scope of application

These instructions for use apply to the following medical devices:

Single use versions of Phaco Sets • Subsonic • F-Sonic • S-Sonic
• Phaco Accessories Set • MediTipS • MediTipC • MediTipF

2. Scope

Intended use: Removal of human lens during eye surgery

Indication: Visual problems that require removal of lens.

Contraindications: There are no specific contraindications for the medical device, apart from general contraindications related to ophthalmic surgery.

User group: The medical device must only be handled by healthcare professionals and applied to the patient by the surgeon. Use by unqualified personnel could result in patient injury.

Patient group: A patient-related, careful preoperative assessment, with well-founded clinical judgment and risk-benefit assessment, is the responsibility of the surgeon.

Clinical benefit: Removal of human lens during eye surgery through small incisions.

Rx ONLY

3. Notes on safety

Warning

Using a damaged or non-sterile medical device can lead to an infection or endophthalmitis of the patient. Do not use if package is damaged or sterile barrier is harmed.

Warning

The medical device must not be used after its expiry date. Use after the expiration date may result in infection / endophthalmitis of the patient.

Warning

Never reuse the medical device. Reuse and / or reesterilization may compromise the device performance and safety, which could cause infections and serious harm to the patient's health and safety.

Never reprocess and re-use the device.

Warning

Pay attention to not touch the activated phaco tip or the phaco tip while disassembling as injuries could occur.

Warning

The medical device may be handled by healthcare professionals and applied to the patient by the surgeon. Use by unqualified personnel could result in patient injury.

Rx ONLY

A patient-related, careful preoperative assessment, with well-

founded clinical judgment and risk-benefit assessment, is the responsibility of the surgeon.

Caution

The User must make sure, that irrigation exits through both irrigation holes

Caution

Make sure to always combine Mediceal components of the correct size

Caution

The modification on the instruments e.g. bending of the tip is not allowed as this could lead to instrument failure (break)

Caution

The User must ensure that appropriate phaco and fluidic parameters are set to avoid damage to the eye.

Caution

The User must make sure that incision size matches the tip / sleeve combination to prevent damage to the patient.

4. Application information

Caution

See illustrated instructions for use. Read and observe the respective surgical system's instructions for use.

5. Disposal

Warning

The medical device may be contaminated with potentially infectious substances of human origin after use. Dispose of the medical device and packaging after use according to the applicable guidelines of biological hazardous waste.

6. Notifications to manufacturer and authorities

Serious incidents occurring in connection with the device must be reported to the manufacturer and the competent authority in which the user and / or patient is established.

7. Symbols not defined by ISO 15223-1:

	Number of medical device in packaging unit		Open here
	To be used by qualified personnel only		Recycling

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