**ENGLISH**

## Medicel VISCOJECT™

### Single-use Lens Injection System

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Intraocular Lens Injection System

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of about 1.5 to 2.8 mm for the very first time, depending on the VISCOJECT™ model and operating technique used, thanks to the new VISCOJECT™ hydraulics system. With the VISCOJECT™ injection system the majority of all foldable, one-piece acrylic lenses with an optical diameter of 6.5 millimeters can be injected, provided that the safe injection of the lens using the particular VISCOJECT™ injection system was validated in advance on the part of the lens manufacturer according to ISO 11979.

The VISCOJECT™ injection system consists of the sterile VISCOJECT™ single-use injector with silicone tip (Fig. 1) and the sterile VISCOJECT™-BIO or VISCOJECT™ single-use cartridge (Fig. 2).

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elastic solutions work for about 30 seconds to achieve sufficient hydration of the lubricating film. After 30 seconds the lens should be injected immediately. Viscoelastic materials may longer their lubricating properties when exposed to air for a longer period of time.

3. Place the foldable lens on the central hinge of the opened cartridge symmetrically as shown in Fig. 4. Align one-piece lenses in the direction of the cartridge regardless of the one-piece design. Asymmetric insertion of the lens can lead to rotation of the lens in the injector.

4. Press the lens into both guide rails with rounded sterile tweezers and gently move the lens back and forth horizontally in order to ensure that it can move freely (Fig. 5).

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silicone tip straight using sterile tweezers. Continue to push the cushion until the inner spring begins to then position the injector. Slowly inject the lens into the eye and simultaneously withdraw the instrument from the eye. To avoid swelling of the silicone tip during ejection from the cartridge, only push the plunger until the lens has emerged completely, even if the plunger is not yet at the limit.

7. Insert the cartridge into the front end of the injector. Push the cartridge into the very front position (Fig. 8).

8. Carefully push the plunger forward and ensure that the silicone tip correctly enters the loading chamber. If the silicone tip cannot be introduced into the loading chamber, retract the plunger to the starting position and align the

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**medicel** SWISS TECHNOLOGY FOR SURGERY

**medicel** SWISS TECHNOLOGY FOR SURGERY

**PLEASE ENGLISH:**  
Do not resterilize / prepare cartridge and injector. Single use only. Reuse and/or resterilization may compromise device performance, which could cause serious harm to the patient's health and safety.

**OPERATIONAL PROCEDURE**  
The appropriate surgical techniques are the responsibility of the respective surgeon. He or she must assess the appropriateness of the relevant procedure based on his or her education and experience.

**GUARANTEE AND LIABILITY LIMITATION**  
The manufacturer guarantees that this product has been manufactured with the appropriate care and assumes no liability for side effects or resulting damages, losses or costs that may arise as a result of the direct or indirect use of this product. Manufacturer's liability is restricted to the performance of repairs resulting from product defects, which are clearly not the result of incorrect handling or the use of lenses not validated with this injector model.

**ATTENTION: US Federal regulations restrict the sale of this product to medical practitioners and those acting on their behalf.\***

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Manufacturer: MEDICEL AG