


# PHACO CASSETTE SYSTEMS

Single-use



Medicel AG

Tel. +41 71 727 10 50

Dornierstrasse 11

Fax +41 71 727 10 55

CH-9423 Altenrhein

info@medicel.com

SWITZERLAND


www.medicel.com

**medicel**  
SWISS TECHNOLOGY FOR SURGERY

## SINGLE-USE CASSETTE SYSTEMS

Apart from a full range of single-use and reusable IA systems, phaco-accessories including phaco tips and infusion sleeves, Medicel offers also a phaco cassette system. The DoubleX K (XXK) is a cassette system for the Series 20000<sup>®</sup> LEGACY<sup>®</sup> - a system for safe, clean and economical procedures.

### Single-use cassette system (XXK)



- Fully compatible with all Series 20000<sup>®</sup> LEGACY<sup>®</sup> cataract systems from Alcon
- Unique total fluidics system which enables ultra-high vacuum techniques
- Absolutely leak-free system

### SINGLE-USE ULTRA-HIGHVAC XXK CASSETTE SYSTEMS, STERILE

consisting of: 1 XXK ULTRA-HIGHVAC cassette, 1 I/A tubing set, 1 Administration set (with filter), 1 Tray support cover

ITEM CODE	PHACO TIPS	PHACO SLEEVES	TEST CHAMBER	INCISION SIZE
<b>XXK100 (Box of 12)</b>	not included	not included	not included	n/a
<b>XXK101 (Box of 12)</b>	not included	19G sleeve, blue	included	2.8mm
<b>XXK120 (Box of 12)</b>	MICROSONIC 20G, 30°	20G sleeve, beige	included	2.5mm