

SWISS COAXIAL SOFT TIP I/A



medicel
swiss evolution

**SWISS COAXIAL I/A System SOFT TIP,
with silicone tip (single-use / sterile)**
consisting of: Coaxial I/A instrument (box of 10)

Item Code	Size	Incision Size*	Aspiration Cannula	Aspiration Port	Aspiration Head
SMS390S	2IG / 0.8mm	2.2mm	angled 45°	0.4mm	silicone tip


* Incision sizes are indications only based on average experiences.
The individually required incision sizes may vary depending on surgical techniques

CE 0482

medicel.com


A **Halma** company

SWISS COAXIAL SOFT TIP I/A




SOFT AND SENSITIVE

Soft contact point to
fragile tissue and IOL

COMPATIBLE

Compatible with standard
luer connector



PERFECT GRIP

Textured handpiece surface for
comfortable handling /grip

GENTLE CONTACT

Silicone Soft Tip
bends upon contact



“ I am impressed
by the design
and apparent
ease of use of
the Soft Tip.”

Dr. Joseph F. Panarelli
USA

The SWISS COAXIAL SOFT TIP I/A allows gentle contact with the sensitive capsular bag tissue & IOL as it is made of soft silicone and visibly bends upon contact. Its textured handpiece grip surface offers comfortable handling and user control. The SWISS COAXIAL SOFT TIP I/A is compatible with common luer connectors. For ease of use, the handpiece is preassembled with sleeve and it requires no sterilization.