



4K F-SONIC



medicel
swiss evolution

4K F-SONIC



High precision phaco tip

The F-SONIC double-step flared design is a high-precision phaco tip. The high irrigation flow results in efficient cooling of the tip, which reduces thermal stress on the surrounding tissue. Maximum freedom through a wide range of sizes compatible with various phaco platforms available on the market.

4K F-SONIC sets (single-use / sterile) consisting of: 10 titanium tips with 10 tip-wrenches / 10 infusion sleeves / 10 test chambers

Item Code	Device Name	Size	Available Degree	Incision Size*
PA603019FK	4K F-Sonic 19G/30°	19G	30° (standard), straight	2.8mm
PA603019CFK	4K F-Sonic 19G/30°, curved	19G	30° (standard), curved	2.8mm
PA603020FK	4K F-Sonic 20G/30°	20G	30° (standard), straight	2.5mm
PA603020CFK	4K F-Sonic 20G/30°, curved	20G	30° (standard), curved	2.5mm
PA603021FK	4K F-Sonic 21G/30°	21G	30° (standard), straight	2.2mm
PA603021CFK	4K F-Sonic 21G/30°, curved	21G	30° (standard), curved	2.2mm
PA603023FK	4K F-Sonic 23G/30°	23G	30° (standard), straight	1.8mm
PA603023CFK	4K F-Sonic 23G/30°, curved	23G	30° (standard), curved	1.8mm

* Incision sizes are indications only based on average experiences. The individually required incision sizes may vary depending on surgical techniques



instruction for use



Medical AG, Dornierstrasse 11, 9423 Altenrhein, +41 71 727 10 50, [medical.com](https://www.medical.com)



digital brochure

made in switzerland