



4K S-SONIC



M
medicel
swiss evolution

4K S-SONIC



Proven S-SONIC phaco tip

S-SONIC is a trusted classic and is available in combination with the perfectly matching silicone sleeves. The high irrigation flow results in efficient cooling of the phacotip, which reduces thermal stress on the surrounding tissue. Maximum freedom through a wide range of sizes compatible with various phaco platforms available on the market.

4K S-SONIC sets (single-use / sterile) consisting of: 10 titanium tips with 10 tip-wrenches / 10 infusion sleeves / 10 test chambers

Item Code	Device Name	Size	Available Degree	Incision Size*
PA603500	4K S-Sonic 19G/30°	19G	30° (standard), straight	2.8mm
PA603500C	4K S-Sonic 19G/30°, curved	19G	30° (standard), curved	2.8mm
PA603510	4K S-Sonic 20G/30°	20G	30° (standard), straight	2.5mm
PA603510C	4K S-Sonic 20G/30°, curved	20G	30° (standard), curved	2.5mm
PA603540	4K S-Sonic 21G/30°	21G	30° (standard), straight	2.2mm
PA603540C	4K S-Sonic 21G/30°, curved	21G	30° (standard), curved	2.2mm

* Incision sizes are indications only based on average experiences. The individually required incision sizes may vary depending on surgical techniques



instruction for use



Medicel AG, Dornierstrasse 11, 9423 Altenrhein, +41 71 727 10 50, [medicel.com](https://www.medicel.com)



digital brochure

made in switzerland