



# ACCUJECT™ SCREW



medicel  
swiss evolution


# ACCUJECT™ SCREW



## Accuject for screw-applications


Accuject screw injectors provides a large range of applications when injecting IOL's. Its versatile design allows for the use of any IOL material and design without limitations on the implant procedure, incision size or the reliable screw application.

### Accuject screw (top-loaded / single-use / sterile)



Item Code	Device Name	Incision size (mm)*			IOL type
		wound-assisted	into-the-wound	into-the-bag	
AST3000	3.0-TL Injector	2.4	2.7 - 2.8	3.1	1-piece lenses
AST2400	2.4-TL Injector	2.1	2.4	2.7	
AST2200	2.2-TL Injector	1.9 - 2.0	2.2	2.5	
AST2100	2.1-TL Injector	1.9	2.1 - 2.2	2.4	
AST2000	2.0-TL Injector	1.8 - 1.9	2.1	2.3 - 2.4	
AST1800	1.8-TL Injector	1.8	2.0	2.2 - 2.3	
AST1600	1.6-TL Injector	1.6	1.8	2.0	hydrophilic only

### Accuject screw (back-loaded / single-use / sterile)



Item Code	Device Name	Incision size (mm)*			IOL type
		wound-assisted	into-the-wound	into-the-bag	
ASB3000	3.0-BL Injector	2.4	2.7 - 2.8	3.1	1-piece lenses
ASB2400	2.4-BL Injector	2.1	2.4	2.7	
ASB2200	2.2-BL Injector	1.9 - 2.0	2.2	2.5	
ASB2100	2.1-BL Injector	1.9	2.1 - 2.2	2.4	
ASB2000	2.0-BL Injector	1.8 - 1.9	2.1	2.3 - 2.4	

\* Incision sizes are indications only based on average experiences. The individually required incision sizes may vary depending on surgical techniques



instruction for use



medicel

Medicel AG, Dornierstrasse 11, 9423 Altenrhein, +41 71 727 10 50, medicel.com



digital brochure

made in switzerland