



The new sustainable bimanual I/A

BICEL terralution combines the best properties of the Medical product ranges currently available in one sustainable product. The track-proven aeroplane shape of the cannulas and the grip structure of the new handpiece design ensure intuitive and comfortable handling as well as precise control during use. It is compatible with all common luer connectors and available in the preferred sizes 23G and 21G.

BICEL terralution I/A (single-use / sterile) consisting of: Bimanual I/A instrument (box of 10)

Item Code	Size	Incision Size*	Aspiration Canula	Aspiration Port	Aspiration Head	Irrigation Canula	Irrigation Port
BIT21A	21G / 0.8mm	1.1mm	aeroplane shape	0.35mm	textured	aeroplane shape	2 ports, 0.45mm
BIT23A	23G / 0.65 mm	0.8mm	aeroplane shape	0.3mm	textured	aeroplane shape	2 ports, 0.35mm

* Incision sizes are indications only based on average experiences. The individually required incision sizes may vary depending on surgical techniques



aeroplane shape