



# ERGOJECT™




**M**  
medicel  
swiss evolution

# ERGOJECT™




## The one-handed screw-injector

Ergoject - the first screw injector for one-handed application. The pen-held design combines the best of both worlds on one-handed (push) and screw. The precise, force-reducing wheel-unit enables controlled delivery and reliable, safe lens injection.



### Ergoject (top-loaded / single-use / sterile)

Item Code	Device Name	Incision size (mm)*			IOL type
		wound-assisted	into-the-wound	into-the-bag	
ERT3000	3.0-TL Injector	2.4	2.7-2.8	3.1	1-piece lenses
ERT2400	2.4-TL Injector	2.1	2.4	2.7	
ERT2200	2.2-TL Injector	1.9-2.0	2.2	2.5	
ERT2100	2.1-TL Injector	1.9	2.1-2.2	2.4	
ERT2000	2.0-TL Injector	1.8-1.9	2.1	2.3-2.4	
ERT1800	1.8-TL Injector	1.8	2.0	2.2-2.3	
ERT1600	1.6-TL Injector	1.6	1.8	2.0	hydrophilic only



### Ergoject (back-loaded / single-use / sterile)

Item Code	Device Name	Incision size (mm)*			IOL type
		wound-assisted	into-the-wound	into-the-bag	
ERB3000	3.0-BL Injector	2.4	2.7 - 2.8	3.1	1-piece lenses
ERB2400	2.4-BL Injector	2.1	2.4	2.7	
ERB2200	2.2-BL Injector	1.9 - 2.0	2.2	2.5	
ERB2100	2.1-BL Injector	1.9	2.1 - 2.2	2.4	
ERB2000	2.0-BL Injector	1.8 - 1.9	2.1	2.3 - 2.4	
ERB1800	1.8-BL Injector	1.8	2.0	2.2-2.3	

\* Incision sizes are indications only based on average experiences. The individually required incision sizes may vary depending on surgical techniques



instruction for use



medical

Medical AG, Dornierstrasse 11, 9423 Altenrhein, +41 71 727 10 50, medical.com



digital brochure

made in switzerland