

VISCOJECT™ -BIO



medicel
swiss evolution

VISCOJECT™ -BIO



A true classic

Provides good injection feel and lens gliding properties.

Viscoject-BIO injector sets (single-use / sterile)

Item Code	Device Name	Incision size (mm)*			IOL type
		wound-assisted	into-the-wound	into-the-bag	
LP604335C	2.7 Injector Set	2.2-2.3	2.5-2.6	2.9	1-piece lenses
LP604340C	2.2 Injector Set	1.9-2.0	2.2	2.5	
LP604350C	1.8 Injector Set	1.8	2.0	2.2	hydrophilic only
LP604370C	1.6 Injector Set	1.6	1.8	2.0	
LP604360C	1.5 Injector Set	1.4	1.6	1.8	for micro-incision lenses only

Viscoject-BIO cartridge sets (single-use / sterile)

Item Code	Device Name	Incision size (mm)*			IOL type
		wound-assisted	into-the-wound	into-the-bag	
LP604235C	2.7 Cartridge Set	2.2-2.3	2.5-2.6	2.9	1-piece lenses
LP604240C	2.2 Cartridge Set	1.9-2.0	2.2	2.5	
LP604250C	1.8 Cartridge Set	1.8	2.0	2.2	hydrophilic only
LP604270C	1.6 Cartridge Set	1.6	1.8	2.0	
LP604260C	1.5 Cartridge Set	1.4	1.6	1.8	

* Incision sizes are indications only based on average experiences. The individually required incision sizes may vary depending on surgical techniques



instruction for use



medicel

Medicel AG, Dornierstrasse 11, 9423 Altenrhein, +41 71 727 10 50, medicel.com



digital brochure

made in switzerland