


VISCOJECT™ ECO



medicel
swiss evolution

VISCOJECT™ ECO



The economic injector

Viscoject eco (single-use / sterile)

Item Code	Device Name	Incision size (mm)*			IOL type
		wound-assisted	into-the-wound	into-the-bag	
VE2700	2.7 Injector Set	2.3	2.6-2.7	2.9-3.0	1-piece lenses
VE2200	2.2 Injector Set	2.1-2.2	2.4	2.7-2.8	

* Incision sizes are indications only based on average experiences. The individually required incision sizes may vary depending on surgical techniques



instruction for use



medicel

Medicel AG, Dornierstrasse 11, 9423 Altenrhein, +41 71 727 10 50, medicel.com



digital brochure

made in switzerland