

# VISCOJECT™ EASY



medicel  
swiss evolution

# VISCOJECT™ EASY



Easy level loading

## Viscoject easy injector sets (single-use / sterile)

Item Code	Device Name	Incision size (mm)*			IOL type
		wound-assisted	into-the-wound	into-the-bag	
VY2900	2.9-IP Injector Set	2.4	2.7	3.0-3.1	1-piece lenses
VY2700	2.7-IP Injector Set	2.2-2.3	2.5-2.6	2.9	
VY2400	2.4-IP Injector Set	2.1	2.4	2.7	
VY2200	2.2-IP Injector Set	1.9-2.0	2.2	2.5	

## Viscoject easy cartridge sets (single-use / sterile)

Item Code	Device Name	Incision size (mm)*			IOL type
		wound-assisted	into-the-wound	into-the-bag	
VC2900	2.9-IP Cartridge Set	2.4	2.7	3.0-3.1	1-piece lenses
VC2700	2.7-IP Cartridge Set	2.2-2.3	2.5-2.6	2.9	
VC2400	2.4-IP Cartridge Set	2.1	2.4	2.7	
VC2200	2.2-IP Cartridge Set	1.9-2.0	2.2	2.5	

\* Incision sizes are indications only based on average experiences. The individually required incision sizes may vary depending on surgical techniques



instruction for use



medicel

Medicel AG, Dornierstrasse 11, 9423 Altenrhein, +41 71 727 10 50, medicel.com



digital brochure

made in switzerland