

ESCRS - Paris 14-18 September 2019 Booth No. E108

The world's leader in IOL injectors and ophthalmic instruments, Medicel, is inviting you to meet with its experts during this year's 37th ESCRS Congress in Paris. Medicel's team will be showcasing the company's revolutionary injector technology, designed to enhance cataract surgery.

Learn more about:

ACCUJECT DUAL™

The world's first 100% 2-in-1 lens injector

The new injector enables surgeons to either operate in push or screw mode, allowing more flexibility. The newest ACCUJECT addition is compatible with all previous ACCUJECT preloaded and non-preloaded systems allowing users of the current ACCUJECT injector to immediately switch to the new injector.



- The ACCUJECT DUAL injector is able to switch from either 100% push to 100% screw at any time of the IOL-injection procedure - without compromise
- Movable finger flanges entirely eliminate usability issues
- Ergonomic design is supported by a unique structured surface which ensures that the ACCUJECT DUAL won't become slippery even when operating it with wet gloves
- Window in the tip of the injector enables the surgeon to have complete visual control of the IOL's travel

**To book an appointment during the exhibition,
please contact Mireille Weder at: m.weder@medicel.com**