

vision



MULTIJECT™

INJECTION SYSTEMS



The reusable solution

Medicel AG

Dornierstrasse 11

CH-9423 Altenrhein

SWITZERLAND

Tel. +41 71 727 10 50

Fax +41 71 727 10 55

info@medicel.com

www.medicel.com

medicel
SWISS TECHNOLOGY FOR SURGERY

MULTIJECT™

The combination of the MULTIJECT™ reusable titanium lens injector and the MULTIJECT™ series cartridges results in a sound and well proven complete lens injection system compatible with silicone as well as acrylic lenses of any design. The MULTIJECT™ allow injection of any style of foldable 1- and 3-piece acrylic and silicone lenses.

Reusable injector



- High grade titanium
- Available as syringe-type or screw-type injector
- To be used with fully single-use cartridge sets
- Hard-tip solution (plunger) for 3-piece IOLs
- Soft-tip solution (with silicone tip) for 1-piece IOLs

Single-use cartridge sets



- Applicable for hydrophilic and hydrophobic lenses
- For 1-piece and 3-piece lenses
- Incision sizes range from 3.2mm to 3.0mm

ITEM CODE	PRODUCT NAME	TYPE	IOL TYPE	TO BE USED WITH MULTIJECT™ CARTRIDGES
LI604200	MULTIJECT™ Injector „GROOVE“	syringe type	1- and 3-piece lenses	LC604215 LC604210 LC604220M
LI604205	MULTIJECT™ Injector „ATRAUMATIC“	syringe type	1-piece lenses	
LI604210	MULTIJECT™ Injector „GROOVE“	screw type	1- and 3-piece lenses	
LI604215	MULTIJECT™ Injector „ATRAUMATIC“	screw type	1-piece lenses	

ITEM CODE	PRODUCT NAME	IOL TYPE	INCISION SIZE
LC604215	MACROGLIDE™ Cartridge	1- and 3-piece lenses	sub 3.5mm
LC604210	OVALGLIDE™ Cartridge	1- and 3-piece lenses	sub 3.2mm
LC604220M	MICROGLIDE™ Cartridge	1- and 3-piece lenses	sub 3.0mm
LS604500	SILICONE CUSHION	to be used with reusable injector LI604205 / LI604215	